



## Portretinis



## Albuminis



## Kvadratas

Lentelė. Apipjovimų dydžiai [Vidus (VID), Išorė (IS), Viršus(VIR), Apačia(AP)]

Puslapis (orientacija)	Įrišimas	Apipjovimas (mm)		
		Tikslus	Standartinis	Presuojant
Portretinis	Kairinis	VID = 2 IS = 2 VIR = 1 AP = 1	VID = 3 IS = 3 VIR = 2 AP = 2	VID = 2 IS = 1 VIR = 1 AP = 1
	Viršutinis	VID = 2 IS = 2 VIR = 1 AP = 1	VID = 3 IS = 3 VIR = 2 AP = 2	VID = 2 IS = 1 VIR = 1 AP = 1
Albuminis	Kairinis	VID = 2 IS = 2 VIR = 1 AP = 1	VID = 4 IS = 4 VIR = 3 AP = 3	VID = 3 IS = 1 VIR = 1 AP = 1
	Viršutinis	VID = 2 IS = 2 VIR = 1 AP = 1	VID = 4 IS = 4 VIR = 3 AP = 3	VID = 3 IS = 1 VIR = 1 AP = 1
Kvadratas	Kairinis	VID = 3 IS = 2 VIR = 2 AP = 2	VID = 4 IS = 3 VIR = 3 AP = 3	VID = 4 IS = 2 VIR = 2 AP = 2
	Viršutinis	VID = 3 IS = 2 VIR = 2 AP = 2	VID = 4 IS = 3 VIR = 3 AP = 3	VID = 4 IS = 2 VIR = 2 AP = 2

*Pastaba.* Spaudos lapo plaušelių kryptis lygiagreti įrišimui.